


TD2823E

Kit Biomateriales Odontológicos

Universidad Andrés Bello campus Viña del Mar

REF		ISO	Grano / Filos	 _{opt.}	 _{max.}
 801.314.012		806.314.001524.012	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 801.314.014		806.314.001524.014	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 830.314.012		806.314.233524.012	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 845KR.314.018		806.314.544524.018	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 846KR.314.016		806.314.545524.016	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 848KR.314.016		806.314.553524.016	Medio (107 µm *)	160.000 min ⁻¹	300.000 min ⁻¹
 8845KR.314.018		806.314.544514.018	Fino (46 µm)	20.000 min ⁻¹	450.000 min ⁻¹
 8862.314.012		806.314.249514.012	Fino (46 µm)	20.000 min ⁻¹	300.000 min ⁻¹
 8368.314.021		806.314.257514.021	Fino (46 µm)	20.000 min ⁻¹	300.000 min ⁻¹
 H1.204.018		500.204.001001.018	-	1.500 min ⁻¹	100.000 min ⁻¹
 9645.204.060		-	-	1.500 min ⁻¹	10.000 min ⁻¹
 310.204		330.204.608000	-	-	30.000 min ⁻¹
 9606.204.060		658.204.030513.060	Medio	6.000 min ⁻¹	15.000 min ⁻¹
 9616.204.060		658.204.030503.060	Fino	6.000 min ⁻¹	15.000 min ⁻¹
 DT-DCP14m		803.204.543523.140	Medio	8.000 min ⁻¹	20.000 min ⁻¹
 DT-DCP14f		803.204.543503.140	Fino	8.000 min ⁻¹	20.000 min ⁻¹
 Enhance Disco		-	Medio	12.000 min ⁻¹	20.000 min ⁻¹

* En algunos instrumentos el tamaño del grano puede variar en función de su forma y tamaño.