

TD2 150D
Kit Prótesis Fija (Preclínico)
Universidad Andrés Bello campus Viña del Mar

REF		ISO	Grano / Filos	 _{opt.}	 _{max.}
  	 6878K.314.016	806.314.298534.016	Grueso (151 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
  	 6879K.314.018	806.314.299534.018	Grueso (151 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
  	 6879K.314.021	806.314.299534.021	Grueso (151 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
  	 8878K.314.016	806.314.298514.016	Fino (46 µm)	20.000 min ⁻¹	450.000 min ⁻¹
  	 8879K.314.018	806.314.299514.018	Fino (46 µm)	20.000 min ⁻¹	450.000 min ⁻¹
  	 8879K.314.021	806.314.299514.021	Fino (46 µm)	20.000 min ⁻¹	450.000 min ⁻¹
 	801.314.018	806.314.001524.018	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
  	 6850.314.012	806.314.199534.012	Grueso (151 µm *)	160.000 min ⁻¹	300.000 min ⁻¹
  	848KR.314.016	806.314.553524.016	Medio (107 µm *)	160.000 min ⁻¹	300.000 min ⁻¹
 	368.314.023	806.314.257524.023	Medio (107 µm *)	160.000 min ⁻¹	300.000 min ⁻¹
 	 8368.314.023	806.314.257514.023	Fino (46 µm)	20.000 min ⁻¹	300.000 min ⁻¹
 	811.314.033	806.314.038524.033	Medio (107 µm *)	100.000 min ⁻¹	100.000 min ⁻¹

* En algunos instrumentos el tamaño del grano puede variar en función de su forma y tamaño.