

TD2200E

Kit de Incrustación (Preclínico)

Universidad Andrés Bello campus Viña del Mar

REF		ISO	Grano / Filos	 opt.	 max.
 845KR.314.016		806.314.544524.016	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 845KR.314.021		806.314.544524.021	Medio (107 µm *)	160.000 min ⁻¹	300.000 min ⁻¹
 846KR.314.016		806.314.545524.016	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 850.314.012		806.314.199524.012	Medio (107 µm *)	160.000 min ⁻¹	300.000 min ⁻¹
 8845KR.314.016		806.314.544514.016	Fino (46 µm)	20.000 min ⁻¹	450.000 min ⁻¹
 8845KR.314.021		806.314.544514.021	Fino (46 µm)	20.000 min ⁻¹	300.000 min ⁻¹
 8846KR.314.016		806.314.545514.016	Fino (46 µm)	20.000 min ⁻¹	450.000 min ⁻¹
 8889.314.009		806.314.540514.009	Fino (46 µm)	20.000 min ⁻¹	300.000 min ⁻¹
 845KRD.314.025		-	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 959KRD.314.018		-	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 8845KR.314.025		806.314.544514.025	Fino (46 µm)	20.000 min ⁻¹	160.000 min ⁻¹
 8959KR.314.018		806.314.584514.018	Fino (46 µm)	20.000 min ⁻¹	450.000 min ⁻¹

* En algunos instrumentos el tamaño del grano puede variar en función de su forma y tamaño.