



TD2199A

Kit Inlays Estéticos Prof. Dr. A. Báez
Universidad Andrés Bello campus Viña del Mar

REF		ISO	Grano / Filos	 _{opt.}	 _{max.}
 845KR.314.018		806.314.544524.018	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 8845KR.314.018		806.314.544514.018	Fino (46 µm)	20.000 min ⁻¹	450.000 min ⁻¹
 845KREF.314.018		806.314.544504.018	Extra fino (25 µm)	20.000 min ⁻¹	450.000 min ⁻¹
 845KRD.314.025		-	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 8845KR.314.025		806.314.544514.025	Fino (46 µm)	20.000 min ⁻¹	160.000 min ⁻¹
 845KREF.314.025		806.314.544504.025	Extra fino (25 µm)	20.000 min ⁻¹	160.000 min ⁻¹
 847KR.314.023		806.314.546524.023	Medio (107 µm *)	160.000 min ⁻¹	300.000 min ⁻¹
 8847KR.314.023		806.314.546514.023	Fino (46 µm)	20.000 min ⁻¹	300.000 min ⁻¹
 10839.314.014		-	Medio (107 µm *)	20.000 min ⁻¹	450.000 min ⁻¹
 959KRD.314.018		-	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 8959KR.314.018		806.314.584514.018	Fino (46 µm)	20.000 min ⁻¹	450.000 min ⁻¹
 8379.314.023		806.314.277514.023	Fino (46 µm)	20.000 min ⁻¹	300.000 min ⁻¹

* En algunos instrumentos el tamaño del grano puede variar en función de su forma y tamaño.