


TD3001C
Kit Carillas Anteriores Prof. Dr. A. Báez
Universidad Andrés Bello campus Viña del Mar

REF		ISO	Grano / Filos	opt. 	max. 
 801.314.014		806.314.001524.014	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 6850.314.014		806.314.199534.014	Grueso (151 µm *)	160.000 min ⁻¹	300.000 min ⁻¹
 6850.314.016		806.314.199534.016	Grueso (151 µm *)	160.000 min ⁻¹	300.000 min ⁻¹
 8850.314.016		806.314.199514.016	Fino (46 µm)	20.000 min ⁻¹	300.000 min ⁻¹
 6850.314.018		806.314.199534.018	Grueso (151 µm *)	160.000 min ⁻¹	300.000 min ⁻¹
 6856.314.021		806.314.198534.021	Grueso (151 µm *)	160.000 min ⁻¹	300.000 min ⁻¹
 8856.314.021		806.314.198514.021	Fino (46 µm)	20.000 min ⁻¹	300.000 min ⁻¹
 868.314.012		806.314.223524.012	Medio (107 µm *)	160.000 min ⁻¹	300.000 min ⁻¹
 8868.314.012		806.314.223514.012	Fino (46 µm)	20.000 min ⁻¹	300.000 min ⁻¹
 8368.314.021		806.314.257514.021	Fino (46 µm)	20.000 min ⁻¹	300.000 min ⁻¹
 H1.204.018		500.204.001001.018	-	1.500 min ⁻¹	100.000 min ⁻¹
 H48L.204.012		500.204.249072.012	12 Filos	20.000 min ⁻¹	300.000 min ⁻¹
 8934A.900.180		-	Fino (46 µm)	5.000 min ⁻¹	25.000 min ⁻¹
 303.204		330.204.603391.050	-	-	15.000 min ⁻¹

* En algunos instrumentos el tamaño del grano puede variar en función de su forma y tamaño.