



TD2245F
Kit Operatoria (Preclínico)
Universidad Andrés Bello campus Viña del Mar

REF	ISO	Grano / Filos	 opt.	 max.
 801.314.010	 806.314.001524.010	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 801.314.012	 806.314.001524.012	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 801.314.014	 806.314.001524.014	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 830.314.010	 806.314.233524.010	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 830.314.014	 806.314.233524.014	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 830.314.016	 806.314.233524.016	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 830L.314.010	 806.314.234524.010	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 830L.314.012	 806.314.234524.012	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 835KR.314.009	 806.314.156524.009	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 845KR.314.016	 806.314.544524.016	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 850.314.014	 806.314.199524.014	Medio (107 µm *)	160.000 min ⁻¹	300.000 min ⁻¹
 8368.314.021	 806.314.257514.021	Fino (46 µm)	20.000 min ⁻¹	300.000 min ⁻¹
 8860.314.012	 806.314.245514.012	Fino (46 µm)	20.000 min ⁻¹	450.000 min ⁻¹
 8863.314.014	 806.314.250514.014	Fino (46 µm)	20.000 min ⁻¹	300.000 min ⁻¹
 8858.314.010	 806.314.165514.010	Fino (46 µm)	20.000 min ⁻¹	300.000 min ⁻¹
 41.204.018	 310.204.001071.018	-	-	40.000 min ⁻¹
 48.204.014	 310.204.243071.014	-	-	40.000 min ⁻¹
 310.204	 330.204.608000	-	-	30.000 min ⁻¹
 9645.204.060	 -	-	1.500 min ⁻¹	10.000 min ⁻¹
 9684.204.040	 -	-	6.000 min ⁻¹	10.000 min ⁻¹

* En algunos instrumentos el tamaño del grano puede variar en función de su forma y tamaño.

TD2245F
Kit Operatoria (Preclínico)
Universidad Andrés Bello campus Viña del Mar

REF	ISO	Grano / Filos	 opt.	 max.
 DT-DCP14m	803.204.543523.140	Medio	8.000 min ⁻¹	20.000 min ⁻¹
 DT-DCP14f	803.204.543503.140	Fino	8.000 min ⁻¹	20.000 min ⁻¹