

TD2135C
Kit Operatoria
Universidad Andrés Bello campus Viña del Mar

REF		ISO	Grano / Filos	 opt.	 max.
 801.314.010		806.314.001524.010	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 801.314.014		806.314.001524.014	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 801.314.018		806.314.001524.018	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 830.314.012		806.314.233524.012	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 830L.314.012		806.314.234524.012	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 838.314.010		806.314.137524.010	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 957AM.314.007		-	Medio (107 µm *)	20.000 min ⁻¹	160.000 min ⁻¹
 H7SM.314.009		-	-	4.000 min ⁻¹	160.000 min ⁻¹
 836KR.314.014		806.314.157524.014	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 880.314.012		806.314.140524.012	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 845KRD.314.025		-	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 8845KR.314.025		806.314.544514.025	Fino (46 µm)	20.000 min ⁻¹	160.000 min ⁻¹
 959KRD.314.018		-	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 8959KR.314.018		806.314.584514.018	Fino (46 µm)	20.000 min ⁻¹	450.000 min ⁻¹
 847KR.314.023		806.314.546524.023	Medio (107 µm *)	160.000 min ⁻¹	300.000 min ⁻¹
 8847KR.314.023		806.314.546514.023	Fino (46 µm)	20.000 min ⁻¹	300.000 min ⁻¹

* En algunos instrumentos el tamaño del grano puede variar en función de su forma y tamaño.