

TD3114A
Kit Endodoncia
Universidad Andrés Bello campus Viña del Mar

REF		ISO	Grano / Filos	 opt.	 max.
 801.314.014		806.314.001524.014	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 801.314.016		806.314.001524.016	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 801.314.018		806.314.001524.018	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 877K.314.014		806.314.297524.014	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 H269QGK.314.016		500.314.219296.016	-	20.000 min ⁻¹	160.000 min ⁻¹
 H1.204.012		500.204.001001.012	-	1.500 min ⁻¹	100.000 min ⁻¹
 H1.204.014		500.204.001001.014	-	1.500 min ⁻¹	100.000 min ⁻¹
 H1.204.016		500.204.001001.016	-	1.500 min ⁻¹	100.000 min ⁻¹

* En algunos instrumentos el tamaño del grano puede variar en función de su forma y tamaño.