

4665.314

Kit Onlay Oclusal según Prof. Dr. Ahlers, Prof. Edelhoff
Universidad Andrés Bello campus Viña del Mar

REF		ISO	Grano / Filos	opt. 	max. 
 855D.314.016	  	-	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 370.314.030		-	Medio (107 µm *)	100.000 min ⁻¹	100.000 min ⁻¹
 370.314.035		-	Medio (107 µm *)	100.000 min ⁻¹	100.000 min ⁻¹
 8849P.314.016	   	-	Fino (46 µm)	20.000 min ⁻¹	160.000 min ⁻¹
 858.314.010	  	806.314.165524.010	Medio (107 µm *)	160.000 min ⁻¹	300.000 min ⁻¹
 8370.314.030		-	Fino (46 µm)	20.000 min ⁻¹	100.000 min ⁻¹
 8370.314.035		-	Fino (46 µm)	20.000 min ⁻¹	100.000 min ⁻¹
 8856.314.014	 	806.314.198514.014	Fino (46 µm)	20.000 min ⁻¹	450.000 min ⁻¹
 8858.314.010	  	806.314.165514.010	Fino (46 µm)	20.000 min ⁻¹	300.000 min ⁻¹

* En algunos instrumentos el tamaño del grano puede variar en función de su forma y tamaño.