

TD2299B
Kit Operatoria Prof. Dr. A. Báez
Universidad Andrés Bello campus Viña del Mar

REF		ISO	Grano / Filos	 opt.	 max.
 801.314.010		806.314.001524.010	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 801.314.014		806.314.001524.014	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 801.314.016		806.314.001524.016	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 830L.314.012		806.314.234524.012	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 838.314.010		806.314.137524.010	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹
 ●●H34L.314.012		500.314.139008.012	-	160.000 min ⁻¹	450.000 min ⁻¹
 ● H35L.314.012		-	-	160.000 min ⁻¹	450.000 min ⁻¹
 ● H856G.314.016		-	-	160.000 min ⁻¹	300.000 min ⁻¹
 801.315.016		806.315.001524.016	Medio (107 µm *)	160.000 min ⁻¹	450.000 min ⁻¹

* En algunos instrumentos el tamaño del grano puede variar en función de su forma y tamaño.