


TD2411G
Kit Pulido Prof. Dr. A. Báez
Universidad Andrés Bello campus Viña del Mar

REF		ISO	Grano / Filos	 opt.	 max.
 8862.314.012		806.314.249514.012	Fino (46 µm)	20.000 min ⁻¹	300.000 min ⁻¹
 8862.314.016		806.314.249514.016	Fino (46 µm)	20.000 min ⁻¹	450.000 min ⁻¹
 8863.314.012		806.314.250514.012	Fino (46 µm)	20.000 min ⁻¹	300.000 min ⁻¹
 6863.314.014		806.314.250534.014	Grueso (151 µm *)	160.000 min ⁻¹	300.000 min ⁻¹
 H48LUF.314.012		500.314.249032.012	30 Filos	20.000 min ⁻¹	300.000 min ⁻¹
 8856.314.016		806.314.198514.016	Fino (46 µm)	20.000 min ⁻¹	450.000 min ⁻¹
 8889.314.009		806.314.540514.009	Fino (46 µm)	20.000 min ⁻¹	300.000 min ⁻¹
 8368.314.023		806.314.257514.023	Fino (46 µm)	20.000 min ⁻¹	300.000 min ⁻¹
 H379Q.314.018		-	-	20.000 min ⁻¹	450.000 min ⁻¹
 H379Q.314.023		-	-	20.000 min ⁻¹	300.000 min ⁻¹
 H41.204.018		500.204.001071.018	12 Filos	20.000 min ⁻¹	100.000 min ⁻¹
 9684.204.040		-	-	6.000 min ⁻¹	10.000 min ⁻¹
 9685.204.060		-	-	6.000 min ⁻¹	10.000 min ⁻¹
 9638.900.190		-	-	6.000 min ⁻¹	10.000 min ⁻¹
 9628.900.220		050.900.373000.220	-	6.000 min ⁻¹	10.000 min ⁻¹
 303.204		330.204.603391.050	-	-	15.000 min ⁻¹
 9645.204.060		-	-	1.500 min ⁻¹	10.000 min ⁻¹
 9645F.204.060		-	-	1.500 min ⁻¹	10.000 min ⁻¹
 310.204		330.204.608000	-	-	30.000 min ⁻¹

* En algunos instrumentos el tamaño del grano puede variar en función de su forma y tamaño.

TD2411G
Kit Pulido Prof. Dr. A. Báez
Universidad Andrés Bello campus Viña del Mar

REF	ISO	Grano / Filos	 opt.	 max.
 DT-DCP14m	 803.204.543523.140	Medio	8.000 min ⁻¹	20.000 min ⁻¹
 DT-DCP14f	 803.204.543503.140	Fino	8.000 min ⁻¹	20.000 min ⁻¹
 Enhance Disco	-	Medio	12.000 min ⁻¹	20.000 min ⁻¹